



For medifa Operating Table Systems

Content

This brochure gives an overview of the complete accessory program for medifa operating table systems. In addition to the complementary accessories for the operating table series MAT5000 BASIC, MAT5000 IV, MOT6000 III and MOT6000S I you will find accessories for nearly each medical application.











| Accessories for MAT5000 I BASIC | 04 | |
|--------------------------------------------|----|--|
| Accessories for MAT5000 IV and MOT6000 III | 09 | |
| Accessories for MOT6000S I | 18 | |
| General accessories | 27 | |
| Head positioning | 37 | |
| Arm positioning | 42 | |
| Leg positioning | 45 | |
| Neurosurgery | 52 | |
| Orthopaedics and traumatology | 54 | |
| Positioning aids | 62 | |

Accessories for MAT5000 I BASIC



Head rest

- Easy adjustment by means of gas spring
- Docking and undocking at back rest
- Radiolucent
- Detachable upholstery

Art.-no. 61148

L x W: 290 x 500 mm

Adjusting range: +45° / -45°

Upholstery: PUR, H = 40 mm

Electrically conductive

Suitable for: 501010, 501020



Horseshoe shaped head rest

- One-piece
- Adjustable by ball joint

Art.-no. 61149

Upholstery: PUR

Electrically conductive

Suitable for: 501010, 501020



Divided leg rest

- Easy adjustment by means of gas spring
- Docking and undocking at pelvis section
- Swivelling
- Radiolucent
- Detachable upholstery
- Pair

Art.-no. 61618

L x W: each 650 x 250 mm


Adjusting range: +20° / -90°

Upholstery: PUR, H = 40 mm

Electrically conductive

Suitable for: 501010, 501020





Leg rest one-piece

- Easy adjustment by means of gas spring
- Docking and undocking at pelvis section
- Radiolucent
- Detachable upholstery

Art.-no. 61628

L x W: 620 x 540 mm

Adjusting range: +20° / -90°

Upholstery: PUR, H = 40 mm

Electrically conductive

Suitable for: 501010, 501020



Carrier

For transport ventilation "Oxylog"

Art.-no. 61091

Attachment at side rail (EU)



Carrier

• For oxygen bottles (2-3 I)

Art.-no. 61090

Accessories for MAT5000 I BASIC



Filing basket

- For attachment at the base of operating table
- Incl. fixation

Art.-no. 61093

Suitable for: 501010, 501020



Floor levelling feet

- Ensures optimal stability of operating table on uneven floors by means of height adjustable levelling feet
- Cannot be retrofitted!

Art.-no. 61094

Suitable for: 501010, 501020



Radiographic cassette carrier

- To insert radiographic cassettes under table top of operating table
- With telescopic rod

Art.-no. 61645

L x W x H: 485 x 410 x 16 mm

Material: Stainless Steel

Suitable for: 501010, 501020





Instrument disposal

Swivelling

Bowl with drainage and tube

Swivelling



Transport box

• For safe transport of operating table

Available upon request!

Art.-no. 61088

L x W x H: 325 x 265 x 65 mm

Material: Stainless steel

Suitable for: 501010, 501020

Art.-no. 61089

L x W x H: 325 x 265 x 65 mm

Material: Stainless steel

Suitable for: 501010, 501020

Art.-no. 61095

L x W x H: 1200 x 800 x 952 mm

Material: Aluminium

Suitable for: 501010, 501020

Accessories for MAT5000 I BASIC



Accessory trolley

- For central storage of accessories
- Five slots for storage of upholstery segements
- Five side rails (25 x 10 mm) for storage of further accessories
- One shelf for storage of small accessories

Art.-no. 341045 1

L x W x H: 765 x 700 x 1670 mm

Four double swivel castors (two with brake)

Suitable for: Upholstery segements of series MAT5000I BASIC as well as accessories for side rails (EU)

Accessories for MAT5000 IV and MOT6000 III





Head rest large

- Easy adjustment by means of gas spring
- Docking and undocking at back rest
- Radiolucent
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61144SF)

Art.-no. 61144

L x W: 285 x 495 mm

Adjusting range: +45° / -45°

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Head rest large

- Easy adjustment by means of gas spring
- Double joint enables additional tilting
- Docking and undocking at back rest
- Radiolucent
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61146SF)

Art.-no. 61146

L x W: 285 x 495 mm

Adjusting range: +45° / -45°

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Head rest small

- Easy adjustment by means of gas spring
- Docking and undocking at back rest
- Radiolucent
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61143SF)

Art.-no. 61143

L x W: 230 x 230 mm

Adjusting range: +45° / -45°

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III

Accessories for MAT5000 IV and MOT6000 III



Upper back rest

- Docking and undocking at lower back rest
- Radiolucent
- Integrated side rails (EU)
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61250SF)

Art.-no. 61250

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: 601700, 601820



Kidney bridge

- Manual adjustable by hand crank
- Cannot be retrofitted!

Art.-no. 61655

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series

MAT5000 IV and 601120



Tripartite back rest

- Easy access to shoulder area by means of three detachable segments
- Docking and undocking at pelvis section (Reverse-Mode)
- Easy adjustment by means of gas spring
- Incl. head rest
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61385SF)

Art.-no. 61385

LxW: each segment 315 x 170 mm

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series

MAT5000 IV and 601120





Tripartite back rest

- Easy access to shoulder area by means of three detachable segements
- Docking and undocking at lower back rest
- Incl. head rest
- Alternatively also available with viscoelastic foam upholsteries, H= 80 mm (Art.-no. 61390SF)


Art.-no. 61390

LxW: each segment 315 x 170 mm

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: 601700, 601820



Pelvis section extension

- Extension of pelvis section by 430 mm
- With cut-out
- Radiolucent
- Integrated side rails (EU)
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61125SF)

Art.-no. 61125

L x W: 430 x 540 mm

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Leg rest

- · Adjustment by means of gas spring
- Docking and undocking at pelvis section
- Radiolucent
- Integrated side rails (EU)
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61625SF)

Art.-no. 61625

L x W: 670 x 540 mm

Adjusting range: +30° / -90°

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III

Accessories for MAT5000 IV and MOT6000 III



Divided leg rest

- · Adjustment by means of gas spring
- Swivelling
- Docking and undocking at pelvis section
- Radiolucent
- Integrated side rails (EU)
- Pair
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61615SF)

Art.-no. 61615

Lx W: each 730 x 250 mm

Adjusting range: +30° / -90°

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Double divided leg rest

- · Adjustment by means of gas spring
- Swivelling and bendable
- Docking and undocking at pelvis section
- Radiolucent
- Integrated side rails (EU)
- Pair
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61660SF)

Art.-no. 61660

L x W: each 720 x 250 mm

Adjusting range: +30° / -90°

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Transfer leg rest

- For short-time positioning of legs
- Docking and undocking at pelvis section

Art.-no. 61635

L x W: 900 x 500 mm

Suitable for: All models of series MAT5000 IV and MOT6000 III





Paediatric surgery rest

- Docking and undocking at pelvis section
- Deatachable upholstery
- Integrated side rails (EU)

Art.-no. 61633

L x W: 700 x 430 mm

Upholstery: Viscoelastic foam, H = 80 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Extension segments

- Segments for extension of operating table's lying surface
- Integrated side rails (EU)
- Pair
- Alternatively also available with viscoelastic foam upholsteries, H= 80 mm (Art.-no. 61355SF)


Art.-no. 61355

LxW: each segment 310 x 105 mm

Upholstery: PUR, H = 60 mm

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Foot switch

- Additional control unit to adjust the main functions height, lateral and Trendelenburg movement
- Tethered
- Cannot be retrofitted!

Art.-no. 604060

Suitable for: All electric models of series MAT5000 IV and all models of series MOT6000 III

Accessories for MAT5000 IV and MOT6000 III



Infrared remote control

- Wireless data transfer
- Control panel for all electric functions of operating table
- Incl. charging station

Art.-no. 604070

Suitable for: All electric models of series MAT5000 IV and all models of series MOT6000 III



Fifth castor

- For easier manoeuvrability and rotatability of operating table
- · Cannot be retrofitted!

Art.-no. 50645

Suitable for: All models of series MAT5000 IV except 505220, 505225, 505320 und 505325



Stabilization cylinder

- Extendible
- For optimal stability of operating table
- Cannot be retrofitted!

Art.-no. 50650

Suitable for: All models of series MOT6000 III





Manual back-up system

- Secures the operability of operating table in case of an electronic defect or breakdown of power supply
- Override of electronic functions by means of manual foot pedal
- Cannot be retrofitted!

Art.-no. 50655

Suitable for: All models of series MOT6000 III



Connecting bar

- Exchange of stainless steel connecting bar at back rest with a radiolucent connecting bar made of special plastic
- Improves radiolucency of back rest

Art.-no. 61670

Suitable for: All models of series MAT5000 IV and MOT6000 III



Universal adapter

 For attachment of a variety of head positioning accessories and neurosurgical accessories

Art.-no. 61100

Suitable for: All models of series MAT5000 IV and MOT6000 III

Accessories for MAT5000 IV and MOT6000 III



Radiographic cassette carrier

- To insert radiographic cassettes under table top of operating table
- With telescopic rod

Art.-no. 61640

L x W x H: 470 x 385 x 16 mm

Material: Stainless steel

Suitable for: All models of series MAT5000 IV and MOT6000 III



Rinsing set

- Incl. bucket, basin, strainer, tube, undercarriage and fixation
- Docking and undocking at pelvis section

Art.-no. 61120

L x W x H (Basin): 325 x 265 x 65 mm

Volume (Bucket): 15 litre

Suitable for: All models of series MAT5000 IV and MOT6000 III



Upholstery

- Alternative covering of lying surface with black viscoelastic foam comfort upholsteries
- With three different foam layers
- With memory effect
- For operating tables with one back section


Art.-no. 5000SF

H = 80 mm

Including:

1 Upholstery for pelvis and back section





Upholstery

- Alternative covering of lying surface with black viscoelastic foam comfort upholsteries
- · With three different foam layers
- · With memory effect
- For operating tables with two back sections

Art.-no. 6000SF

H = 80 mm

Including:

- 1 Upholstery for pelvis section
- 1 Upholstery for lower back section



Accessory trolley

- For central storage of accessories
- Five slots for storage of upholstery segements
- Five side rails (25 x 10 mm) for storage of further accessories
- One shelf for storage of small accessories


Art.-no. 341045 2

L x W x H: 765 x 700 x 1670 mm

Four double swivel castors (two with brake)

Suitable for: Upholstery segements of series MAT5000 IV and MOT6000 III as well as accessories for side rails (EU)

Accessories for MOT6000S I



Head rest

- Easy adjustment by means of gas spring
- Docking and undocking at back rest
- Radiolucent
- Detachable upholstery enables easy cleaning and disinfection

Art.-no. 63144

L x W: 266 x 550 mm

Adjusting range: -45° / +45°


Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603100, 603200,

603300



Head rest

- Easy adjustment by means of gas spring
- Double joint enables additional tilting
- Docking and undocking at back rest
- Radiolucent
- Detachable upholstery enables easy cleaning and disinfection

Art.-no. 63146

L x W: 266 x 550 mm

Adjusting range: -45° / +45°


Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603100, 603200,

603300



Side rails at head rest

- Additional side rail at head-end of head rest
- Piece

Side rails at head rest

- Additional side rails at both sides of head rest
- Pair


Art.-no. 63148

Suitable for: 63144, 63146

Art.-no. 63149

Suitable for: 63144, 63146





Upper back rest

- Docking and undocking at lower back rest
- Radiolucent
- Detachable upholstery enables easy cleaning and disinfection
- Integrated side rails (EU)

Art.-no. 63250 1

L x W: 268 x 550 mm


Adjusting range: -35° / +50

Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603100



Upper back rest

- Docking and undocking at lower back rest
- Radiolucent
- Detachable upholstery enables easy cleaning and disinfection
- Integrated side rails (EU)

Art.-no. 63250 2

L x W: 268 x 550 mm

Adjusting range: -35° / +50

Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603200, 603300



Tripartite back rest

- Easy access to shoulder area by means of three detachable segments
- Docking and undocking at lower back rest
- Radiolucent
- Detachable upholsteries enable easy cleaning and disinfection
- Incl. head rest

Art.-no. 63390

L x W: 320 x 550 mm

Adjusting range: -35° / +50

Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603100, 603200,

603300

Accessories for MOT6000S I



Pelvis section extension

- Extension of pelvis section by 400 mm
- With cut-out
- Radiolucent
- Detachable upholstery enables easy cleaning and disinfection
- Integrated side rails (EU)

Art.-no. 63125_1

L x W: 400 x 550 mm

Upholstery: Viscoelastic foam, H = 80 mm

Electrically conductive

Suitable for: 603100



Pelvis section extension

- Extension of pelvis section by 400 mm
- With cut-out
- Radiolucent
- Detachable upholstery enables easy cleaning and disinfection
- Integrated side rails (EU)


Art.-no. 63125 2

L x W: 400 x 550 mm

Upholstery: Viscoelastic foam, H = 80 mm

Electrically conductive

Suitable for: 603200, 603300



Leg rest

- Electrically adjustable
- Docking and undocking at pelvis section
- Radiolucent
- Detachable upholstery enables easy cleaning and disinfection
- Integrated side rails (EU)

Art.-no. 63625

L x W: 660 x 550 mm

Adjusting range: -90° / +70°


Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603200, 603300





Divided leg rest

- · Easy adjustment by means of gas spring
- Swivelling
- Docking and undocking at pelvis section
- Radiolucent
- Detachable upholsteries enable easy cleaning and disinfection
- Integrated side rails (EU)
- Pair



- Electrically adjustable
- Swivelling
- Docking and undocking at pelvis section
- Radiolucent
- Detachable upholsteries enable easy cleaning and disinfection
- Integrated side rails (EU)
- Pair

Art.-no. 63615 1

L x W: each 660 x 266 mm


Adjusting range: -90° / +30°

Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603100



Art.-no. 63615 2

L x W: each 660 x 266 mm

Adjusting range: -90° / +80°

 $Uphol stery: Viscoelastic \ foam,$

H = 80 mm

Electrically conductive

Suitable for: 603200, 603300



Double divided leg rest

- · Electrically adjustable
- Swivelling and bendable
- Docking and undocking at pelvis section
- Detachable upholsteries enable easy cleaning and disinfection
- Integrated side rails (EU)
- Pair

Art.-no. 63660

L x W: each 662 x 266 mm

Adjusting range: -90° / +80°

Upholstery: Viscoelastic foam,

H = 80 mm

Electrically conductive

Suitable for: 603200, 603300

Accessories for MOT6000S I



Carbon positioning rest

- For intraoperative x-ray
- Optimal radiolucency (360°)
- Docking and undocking at pelvis section or lower back rest
- Incl. upholstery

Item in peperation and currently not available!

Art.-no. 63630

Lx W: 900 x 540 mm

Upholstery: Viscoelastic foam, H= 80 mm

Suitable for: 603100, 603200, 603300



Transfer leg rest

- For short-time positioning of legs
- Docking and undocking at pelvis section

Item in peperation and currently not available!

Art.-no. 63635

L x W: 900 x 500 mm

Suitable for: 603100, 603200, 603300



Extension segments

- For extension of the OT-Table's lying surface
- Integrated side rails (EU)
- Pair

Art.-no. 61355SF

L x W: each 310 x 110 mm

Upholstery: Viscoelastic foam, H = 80 mm

Electrically conductive

Suitable for: 603100, 603200,

603300





Foot switch

- Additional control unit to adjust the main functions height, lateral and Trendelenburg movement
- Tethered
- Cannot be retrofitted!

Art.-no. 634060

Suitable for: 603100, 603200, 603300




Infrared remote control

- Wireless data transfer
- Control panel for all electric functions of operating table
- Incl. charging station

Art.-no. 634070

Suitable for: 603100, 603200, 603300



Sliding

- Electric sliding of table top up to 400 mm
- Cannot be retrofitted!

Art.-no. 63640

Suitable for: 603200, 603300

Accessories for MOT6000S I



Back-up system

- System ensures the operability of operating table's basic functions (height, lateral, Trendelenburg and brakes) in case of electronic defect by means of 24V battery
- Cannot be retrofitted!

Art.-no. 63655

Suitable for: 603100, 603200, 603300



Electronic drive

- Motorized fith castor for manoeuvring and positioning of operating table
- Cannot be retrofitted!

Art.-no. 63646

Suitable for: 603100, 603200, 603300



Hand grips for electronic drive

- For manoeuvring and positioning of operating table by means of electronic drive at operating table's head end
- With integrated speed control
- Movement in forward and backward direction possible

Art.-no. 63570

Attachment at side rail (EU)

Suitable for: 63646





Radiographic cassette tunnel

• To insert radiographic cassette carrier

uitable for 603100

Suitable for: 603100

Radiographic cassette tunnel

• To insert radiographic cassette carrier

Art.-no. 63642 2

Art.-no. 63642 1


Suitable for: 603200

Radiographic cassette tunnel

• To insert radiographic cassette carrier

Art.-no. 63642 3

Suitable for: 603300




Radiographic cassette carrier

- To insert radiographic cassettes under table top of operating table
- With telescopic rod

Art.-no. 63643

L x W x H: 470 x 385 x 16 mm

Material: Stainless steel Suitable for: 63642



Universal adapter

 For attachment of a variety of head positioning accessories and neurosurgical accessories

Art.-no. 63100

Suitable for: 603100, 603200, 603300

Accessories for MOT60005 I



Accessory trolley

- For central storage of accessories
- Five slots for storage of upholstery segements
- Five side rails (25 x 10 mm) for storage of further accessories
- One shelf for storage of small accessories

Art.-no. 341045 3

L x W x H: 765 x 700 x 1670 mm

Four double swivel castors (two with brake)

Suitable for: Upholstery segements of series MOT6000S 1 as well as accessories for side rails (EU)

General Accessories





Clamp

- For fixation of accessories with round rod (18 mm)
- Rotatable

Art.-no. 61511

Attachment at side rail (EU)



Clamp

- For fixation of accessories with round rod (18 mm)
- With quick lock and release lever
- Rotatable

Art.-no. 61<u>512</u>

Attachment at side rail (EU)



Clamp

 For fixation of accessories with square rod (25 x 10 mm)

Art.-no. 61551

General Accessories



Clamp

 For fixation of accessories with square rod (20 x 20 mm)

Art.-no. 61552

Attachment at side rail (EU)



Clamp

 For fixation of accessories with round rod (18 mm)

Art.-no. 61553

Attachment at side rail (EU)



Clamp

- For fixation of accessories with round rod (18 mm)
- With quick lock and release lever

Art.-no. 61554





Anaesthesia screen

- Extendible to both sides
- Incl. clamp

Art.-no. 61750

Attachment at side rail (EU)



Arm restraint cuff

• For attachment at anaesthesia screen

Art.-no. 61755

Suitable for: 61750



Flexible anaesthesia screen

- Bendable
- · Individual adjustable
- Incl. clamp

Art.-no. 61751

L: 2000 mm

General Accessories



Tube support

- For fixation of breathing tube
- Bendable support is infinitely variable
- Adjustable in height
- Incl. clamp

Art.-no. 61752

Attachment at side rail (EU)




Infusion holder

- Extendible
- Four hooks
- Incl. clamp

Art.-no. 61800

Attachment at side rail (EU)



Infusion holder

- Extendible
- Rotatable
- Four hooks
- Incl. rotatable clamp

Art.-no. 61805





Mayo table

- Adjustable in height
- Rotatable
- Inclinable
- Incl. clamp

Art.-no. 61410

L x W: 450 x 350 mm

Attachment at side rail (EU)



Side rail extension

- Extends side rail up to 320 mm
- For attachement of accessories at side rail (EU)
- Incl. clamp

Art.-no. 61560

Attachment at side rail (EU)



Sliding handles

- For easy and comfortable movement of operating table
- Incl. clamp
- Pair

Art.-no. 61571

General Accessories



Side support

- Incl. clamp
- Pair

Art.-no. 61341

L x W: each 220 x 110 mm

Upholstery: PUR

Electrically conductive

Attachment at side rail (EU)



Shoulder support

- Adjustable in height
- Extendible
- Incl. clamp
- Pair

Art.-no. 61351

L x W: each 220 x 110 mm

Upholstery: PUR

Electrically conductive

Attachment at side rail (EU)



Double side support

- Support of pelvis section
- Adjustable in height
- Extendible
- Incl. clamp

Art.-no. 61362

Upholstery: PUR, Ø: 100 mm

Electrically conductive





Holder

- For attachment of body supports
- · Adjustable in height
- Extendible
- Swivelling (360°)
- Positioning of body supports by means of ball joint

Art.-no. 61361

Attachment at side rail (EU)



Joint arm

- For attachment of body supports
- Flexible positioning of body supports by means of three joints
- Quick lock and release of joint arm by means of central screw

Art.-no. 61365

Attachment at side rail (EU)



Body support

To support pubis, sacrum and sternum

Art.-no. 61361 1

Ø: 100 mm

Upholstery: PUR

Electrically conductive

Suitable for: 61361, 61365

General Accessories



Body support

• To support side and shoulder

Art.-no. 61361_2

L x W: 220 x 110 mm

Upholstery: PUR

Electrically conductive

Suitable for: 61361, 61365



Body support

• To support back and buttocks

Art.-no. 61361 3

L x W: 220 x 110 mm

Upholstery: PUR

Electrically conductive

Suitable for: 61361, 61365



Body support

Supporting roll for lateral positioning of patient

Art.-no. 61361_4

L x Ø: 190 x 100 mm

Upholstery: PUR

Electrically conductive

Suitable for: 61361, 61365





Instrument disposal

- Round form
- Swivelling
- Incl. clamp

Art.-no. 91611

 $\emptyset = 280 \, \text{mm}$

Material: Stainless steel

Attachment at side rail (EU)



Instrument disposal

- Squared form
- Swivelling
- Incl. clamp

Art.-no. 91612

L x W x H: 325 x 265 x 65 mm

Material: Stainless steel

Attachment at side rail (EU)



Bowl with drainage and tube

- For the intake and drainage of fluids
- Swivelling
- Incl. clamp

Art.-no. 91613

 $\emptyset = 320 \, \text{mm}$

Material (Bowl): Stainless steel

General Accessories



Bucket

With lid and tube-inlet

Art.-no. 91608

Volumen: 12 Liter

Material: Stainless Steel

Suitable for: 91613



Hand switch support

- To place hand switch at lateral or head-sided side rail
- Tiltable

Art.-no. 61575

Material: Stainless steel

Head positioning






Universal adapter

 For attachment of a variety of head positioning accessories and neurosurgical accessories

Art.-no. 61100

Suitable for: All models of series MAT5000 IV und MOT6000 III




Universal adapter

 For attachment of a variety of head positioning accessories and neurosurgical accessories

Art.-no. 63100

Suitable for: 603100, 603200, 603300



Horseshoe shaped head rest

- For positioning head in prone or supine position
- · Oval form with cut-out
- Tiltable and rotatable by means of ball joint
- · Adjustable in height

Art.-no. 61114_1

Upholstery: PUR

Electrically conductive

Attachment at 61100, 61385, 61390, 63100 or 63390

Head positioning



Horseshoe shaped head rest

- For positioning head in supine position
- Rotatable and tiltable by means of ball joint
- Adjustable in height

Art.-no. 61114 2

Upholstery: PUR

Electrically conductive

Attachment at 61100, 61385, 61390, 63100 or 63390



Horseshoe shaped head rest

- For positioning head in prone or supine position
- Two-part
- Extendible in width
- Rotatable and tiltable by means of ball joint
- Incl. detachable upholsteries and fixture strap

Art.-no. 61114 3

Upholstery: Viscoelastic foam

Electrically conductive

Attachment at 61100, 61385, 61390, 63100 or 63390



Head rest

- Easy adjustment by means of gas spring
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61116SF)

Art.-no. 61116

L x W: 230 x 230 mm

Upholstery: PUR, H = 60 mm:

Electrically conductive

Attachment at 61100, 61385, 61390, 63100 or 63390





Head support

- For comfortable and safe positioning of patient's head during shoulder surgeries
- Attachment at tripartite back rest or universal adapter
- Extendible
- Incl. upholstery

Art.-no. 61395

Upholstery: Viscoelastic foam

Electrically conductive

Attachment at 61100, 61385, 61390, 63100 or 63390



Joint arm

- For attachment of head positioning accessories
- Flexible positioning of body supports by means of three joints
- Quick lock and release of joint arm by means of central screw

Art.-no. 61117

Attachment at 61100 or 63100



Head rest

- For comfortable and safe positioning of patient's head
- Alternatively also available with viscoelastic foam upholstery, H= 80 mm (Art.-no. 61116SF)

Art.-no. 61117 1

L x W: 230 x 230 mm

Upholstery: PUR, H= 60 mm

Electrically conductive

Attachment at 61117

Head positioning



Head support

- For stable positioning and fixation of patient's head
- Incl. upholsteries and fixation strap

Art.-no. 61117 2

Upholstery: Viscoelastic foam Electrically conductive Attachment at 61117



Horseshoe shaped head rest

- For positioning head in prone or supine position
- Oval form
- With cut-out

Art.-no. 61117 3

Upholstery: PUR

Electrically conductive

Attachment at 61117



Horseshoe shaped head rest

For positioning head in supine position

Art.-no. 61117 4

Upholstery: PUR

Electrically conductive

Attachment at 61117





Horseshoe shaped head rest

- For positioning head in prone or supine position
- Two-part
- Extendible in width
- Incl. detachable upholsteries and fixture strap

Art.-no. 61117 5

Upholstery: Viscoelastic foam

Electrically conductive Attachment at 61117



Wrist support ring

With connecting piece for horseshoe shaped head rest

Art.-no. 61111

Material: Stainless steel

Suitable for: 61114_2



Wrist support ring

With connecting piece for horseshoe shaped head rest

Art.-no. 61113

Material: Stainless steel

Suitable for: 61114_1, 61114_3, 61117 3, 61117 4, 61117 5

Arm positioning



Arm rest

- Rotatable and tiltable by means of ball joint
- Adjustable in height
- · Horizontally slidable
- Incl. clamp and fixture straps

Art.-no. 61211

L x W: 400 x 115 mm

Upholstery: PUR, H = 40 mm

Electrically conductive

Attachment at side rail (EU)



Intraoperative arm positioner

- For positioning of patient's arm during shoulder surgeries in beach chair position
- Quick lock and release of joint arm by means of central screw
- Incl. clamp and fixture straps

Art.-no. 61216

L x W: 400 x 115 mm

Upholstery: PUR, H = 40 mm

Electrically conductive

Attachment at side rail (EU)



Arm rest

- Three joints enable flexible positioning of patient's arm
- Quick lock and release of joint arm by means of central screw
- Incl. clamp and fixture straps

Art.-no. 61217

L x W: 400 x 115 mm

Upholstery: PUR, H = 40 mm

Electrically conductive





Arm rest

- Rotatable by 360°, swivelling and tiltable by means of one-hand control
- Adjustable in height
- Incl. clamp and fixture straps

Art.-no. 61221

LxW: 600 x 165 mm

Upholstery: PUR, H = 60 mm

Electrically conductive

Attachment at side rail (EU)



Arm rest

- Rotatable by 360°, swivelling and tiltable by means of one-hand control
- · Adjustable in height
- Detachable upholstery enables easy cleaning and disinfection
- Incl. clamp and fixture straps

Art.-no. 61222

LxW: 600 x 165 mm

Upholstery: Viscoelastic foam, H = 60 mm

Electrically conductive

Attachment at side rail (EU)



Arm rest

- Rotatable by 360°, swivelling and tiltable by means of one-hand control
- Adjustable in height
- Detachable upholstery enables easy cleaning and disinfection
- Incl. clamp and fixture straps

Art.-no. 61223

LxW: 450 x 165 mm

Upholstery: Viscoelastic foam, H = 60 mm

Electrically conductive

Arm positioning



Arm rest

- Swivelling
- Incl. clamp and fixture strap


Art.-no. 61220

LxW: 600 x 165 mm

Upholstery: PUR, H = 60 mm

Electrically conductive

Attachment at side rail (EU)



Side support / arm rest

- For stabilisation of patient or positioning of arm
- Hinged in three steps (upwards, horizontally, downwards)

Art.-no. 61092

L x H: 600 x 140 / 235 mm

Attachment at side rail (EU)



Upholstery

• For side support / arm rest 61092

Art.-no. 56031

Upholstery: Artificial leather covered foam, H= 20 mm

Electrically conductive





Hand cuffs

- For fixation of hands or arms
- Incl. clamp
- Pair

Art.-no. 61400

L x W: each 180 x 80 mm

Leg positioning



Leg rest according to Goepel

- For positioning patient in lithotomy position
- Piece

Art.-no. 61700

Upholstery: PUR

Electrically conductive

Attachment at side rail (EU)

Rotatable clamp (61511 or 61512) and fixture straps (B81600M1) to be ordered separately



Leg rest according to Goepel Comfort

- For positioning patient in lithotomy position
- With grey foamed PUR knee bolsters
- Pair

Art.-no. 81705

Electrically conductive

Attachment at side rail (EU)

Rotatable clamp (61511 or 61512) and fixture straps (B81600M1) to be ordered separately



Leg rest

- For positioning patient in lithotomy position
- Stepless positioning of legs
- Easy adjustment by means of gas spring supported one-hand control
- Incl. upholsteries and clamps
- Pair

Art.-no. 61720

Max. patient weight: 159 kg Attachment at side rail (EU)





Leg rest

- Especially for the positioning of adipose patients in lithotomy position
- · Tiltable by means of ball joint
- Incl. upholstery and fixture strap
- Pair

Art.-no. 61710

Upholstery: 300 x 285 x 155 / 100 mm

Electrically conductive

Attachment at side rail (EU)

Rotatable Clamp (61511 or 61512) to be ordered separately



Elbow support

- To support elbows of surgeon during procedures in lithotomy position
- Incl. clamp
- Pair

Art.-no. 61780

L x W: 145 x 120 mm

Upholstery: PUR

Electrically conductive

Attachment at side rail (EU)



Lower leg fixation

- Incl. clamp
- Piece

Art.-no. 61760

Upholstery: PUR

Electrically conductive

Leg positioning



Thigh fixation

- Incl. clamp
- Piece

Art.-no. 61770

Upholstery: PUR

Electrically conductive

Attachment at side rail (EU)



Rectal positioning device

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- Adjustable in height (200 mm)
- Integrated side rails (EU)

Art.-no. 61632

L x W x H: 370 x 500 x 490 - 690 mm

Max. patient weight: 325 kg

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Rectal positioning device

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- Adjustable in height (200 mm)
- Integrated side rails (EU)

Art.-no. 63632

L x W x H: 370 x 500 x 490 - 690 mm

Max. patient weight: 325 kg

Electrically conductive

Suitable for: All models of series

MOT6000S I





Knee supports

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- · Adjustable in height
- Swivelling
- Piece

Art.-no. 61150

Upholstery: PUR

Electrically conductive

Suitable for: All models of series MAT5000 IV and MOT6000 III



Rectoscopy roll

- Adjustable in height
- Incl. clamp

Art.-no. 61315

Upholstery: PUR

Electrically conductive

Attachment at side rail (EU)



Knee support

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- Incl. clamp

Art.-no. 61630

L x W: 525 x 245 mm

Upholstery: PUR

Electrically conductive

Leg positioning



Foot plates

- To stabilize patient in Anti-Trendelenburg position
- For extension of lying surface (foot plates can also be hinged down)
- Incl. clamps
- Pair

Art.-no. 61631

L x W: each 270 x 270 mm

Upholstery: PUR, H = 60 mm

Electrically conductive

Attachment at side rail (EU)



Body and leg strap

- Steplessly adjustable
- Incl. clamp

Art.-no. 61405

L x W: 2300 x 80 mm

Attachment at side rail (EU)



Leg strap

- For fixation of legs at divided leg rest
- Steplessly adjustable
- Pair

Art.-no. 61406

Suitable for: 61618, 61615, 61660, 63615_1, 63615_2, 63660_1, 63660_2





Leg straps

- For divided leg rests
- Incl. clamp
- Pair

Art.-no. 61407

L: each 1200 mm

Safety load: 454 kg (when using as pair)

Attachment at side rail (EU)

Suitable for: 61618, 61615, 61660, 63615_1, 63615_2, 63660_1, 63660_2

Neurosurgery



DORO® Standard holder

- For fixation of skull clamp or horseshoe headrest
- Adjustable with one hand
- · Large adjustment radius
- · Very low weight

Art.-no. 661065

Only in connection with 61100 or 63100



DORO® Swivel adapter

- Connects standard holder (661065) with skull clamp (661069)
- Enables rotation adjustment of skull clamp

Art.-no. 661067

Suitable for: 661065, 661069, 661076 and 661078



DORO® Skull clamp

- For fixation of head during neurosurgical procedures
- Locking mechanism enables quick and secure fixation

Art.-no. 661069

Only in connection with 661065 and 661067

Skull clamps (661073 and 661074) are not included in delivery!





DORO® Skull pins adults

- · Can be sterilised
- For multiple uses
- Three pieces

DORO® Skull pins children

- Can be sterilised.
- For multiple uses
- Three pieces

DORO® Cross bar adapter

- For secure fixation of head in seated position of patient
- Adaptable to the width of operating table
- Incl. clamps

Art.-no. 661073

Material: Stainless steel Suitable for 661069

Art.-no. 661074

Material: Stainless steel Suitable for 661069

Art.-no. 661071

Suitable for: 661065

Attachment at side rail (EU)



DORO® Horseshoe headrest for adults

- For positioning of head in prone or supine position
- Incl. detachable gel pads

- DORO® Horseshoe headrest for children
- For positioning of head in prone or supine position
- Incl. detachable gel pads

Art.-no. 661076

Only in connection with 661065 and 661067

Art.-no. 661078

Only in connection with 661065 and 661067



Orthopaedics and traumatology

Extension device

- For attachment at operating table
- Levitating
- Max. patient weight:
 - Levitating: 136 kg
 - With supporting bars: 250 kg

Included in delivery:

- Fixation adapter for operating table incl. pelvis part with upholstery
- Countertraction post for perineum and for neck of femur
- Two swivelling extension bars with double joints
- Two swivelling and height adjustable drawing spindles with adaption for extension shoe holders
- One trolley with storage basket for additional accessories
- Two height adjustable supporting bars



Art.-no. 63000 1

Material: All metal parts made of stainless steel

Suitable for: All models of series MAT5000 IV and MOT6000 III

Alternatively also available with viscoelastic foam upholstery, H = 80 mm (Art.-no. 63000 1SF)

Art.-no. 63000 2

Material: All metal parts made of stainless steel

Suitable for: All models of series MOT6000S I







Extension shoe for adults

- · For fixation of foot
- Piece

Art.-no. 61010

Suitable for: 63000 1, 63000 2

Extension shoe for children

- For fixation of foot
- Piece

Art.-no. 61020

Suitable for: 63000 1, 63000 2



Countertraction post for femur

- For extension in lateral position
- Radiolucent perinneal arc

Art.-no. 63025 1

Upholstery: PUR

Suitable for: 63000 1

Countertraction post for femur

- For extension in lateral position
- Radiolucent perinneal arc

Art.-no. 63025 2

Upholstery: PUR

Suitable for: 63000 2



Countertraction post for tibia and fibula

- For positioning and fixation of patient during surgeries at lower leg
- Radiolucent padded roll
- Adjustable in height
- Inclinable

Art.-no. 63045

Material: Metal parts made of

stainless steel

Upholstery: PUR

Suitable for: 63000 1, 63000 2

Orthopaedics and traumatology



Pelvis plaster cast unit

- Usable for adults and children
- Slidable countertraction post

Art.-no. 63030

Material: Metal parts made of stainless steel

Upholstery: PUR

Suitable for: 63000_1, 63000_2



Horseshoe holder "Kirschner"

· For extensions using Steinmann nails

Art.-no. 63060

Material: Stainless steel

Suitable for: 63000_1, 63000_2



Leg support

- For positioning of healthy leg during extension
- Tiltable upholstery by means of ball joint
- Incl. fixture straps

Art.-no. 63026

Material: Metal parts made of stainless steel

Upholstery: PUR

Suitable for: 63000_1, 63000_2





Transfer leg rest

- For preoperative positioning of legs
- Incl. upholsteries (H= 35 mm)
- Pair

Transfer leg rest

- For preoperative positioning of legs
- Incl. upholsteries (H= 35 mm)
- Pair



Art.-no. 63070

L x W: 600 x 200 mm Suitable for: 63000 1

Art.-no. 63070 2

L x W: 600 x 200 mm Suitable for: 63000 2

Trolley for extension device

 For storage and transport of extension device and accessories

Art.-no. 63090

Material: Stainless steel

Suitable for: 63000_1, 63000_2

Orthopaedics and traumatology



Arm extension device

- For intramedullary nailing of humerus
- Radiolucent countertraction post
- Incl. tension rod, countertraction post and clamp

Art.-no. 62050

Attachment at side rail (EU)

Upholstery: PUR




Weinberger hand cuff

- Secure fixation of hand and upper arm during intramedullary nailing of humerus
- Radiolucent
- Can be sterilised

Art.-no. 62055

Suitable for: 62050



Upper arm support

- For positioning of upper arm
- Radiolucent
- Adjustable in height

Art.-no. 61226

L x W: 280 x 280 mm

Upholstery: Artificial leather covered

foam, H = 30 mm

Electrically conductive





Arm and hand surgical table

- · Adjustable in height
- Levitating
- Radiolucent
- Incl. clamps and upholstery

Art.-no. 61231

LxW: 600 x 340 mm

Upholstery: Viscoelastic foam, H = 35 mm

Electrically conductive

Attachment at side rail (EU)



Arm and hand surgical table

- Adjustable in height
- Radiolucent
- With tunnel for radiographic cassette carrier
- Incl. clamp, supporting leg and upholstery

Art.-no. 61232

L x W: 645 x 295 mm

Upholstery: Artificial leather covered foam, H = 25 mm

Electrically conductive

Attachment at side rail (EU)



Arm and hand surgical table

- Table top made of carbon
- Radiolucent
- · Adjustable in height
- Incl. clamp, upholstery and supporting leg

Art.-no. 56230

L x W: 815 x 430 mm

Upholstery: Viscoelastic foam, H = 30 mm

Electrically conductive

Orthopaedics and traumatology



Meniscotomy roll

- Adjustable in height
- Incl. clamp

Art.-no. 61330

Upholstery: PUR

Electrically conductive

Attachment at side rail (EU)



Arthroscopic leg holder system

- For fixation of thigh
- Horizontally and vertically adjustable
- Width of leg holder adjustable in three steps
- Incl.upholstery, clamp and fixture strap

Art.-no. 61328

Upholstery: Artificial leather covered foam

Electrically conductive





Knee positioning device

- Electrically adjustable by means of foot switch
- Intraoperative flexion and extension of knee
- Fixation of thigh
- Width of leg holder adjustable in three steps
- Suitable for different leg lengths
- Incl. upholstery, clamps and fixture strap

Art.-no. 61327

Tilt angle: 0° - 60°

Upholstery: Artificial leather covered

foam

Electrically conductive

Attachment at side rail (EU)

Power supply: 230 V

Battery: 24 V, 1.2 Ah

Positioning aids



Head support

- For positioning of patient in supine position
- Protects from positioning injuries

Art.-no. 61147

L x W: 250 x 250 mm

Material: PUR

Electrically conductive



Head ring

- For positioning of patient in supine position
- Protects from positioning injuries

Art.-no. 61115

H: 45 mm

Ø (inner): 75 mm

Ø (outer): 200 mm

Material: Silicone gel



Head positioning pad adults

- For positioning of patient in prone position
- Protects from positioning injuries

Head positioning pad children

- For positioning of patient in prone position
- Protects from positioning injuries

Art.-no. 61910

L x W x H: 280 x 240 x 150 mm


Material: Silicone gel

Art.-no. 61911

L x W x H: 260 x 230 x 140 mm

Material: Silicone gel





Head positioning pad


- For positioning patient in supine position
- · Protects from positioning injuries

Art.-no. 61905

L x W x H: 400 x 500 x 100 mm

Material: Artificial leather covered foam

Electrically conductive



Lateral positioning pad

- For stable positioning during hip, shoulder or kidney surgeries
- Cut-out to relieve pressure from shoulder joint and upper arm
- Protects from positioning injuries

Art.-no. 61915

L x W x H: 700 x 500 x 160 / 25 mm

Material: PU coated foam

Electrically conductive



Spinal disc pad

- For positioning during surgeries at spinal disc
- Radiolucent

Art.-no. 61900

L x W x H: 820 x 450 x 90 mm

Material: Viscoelastic foam

Electrically conductive

Positioning aids



Flat bottom chest roll large

 Protects from positioning injuries at legs or chest

Art.-no. 61922

 $L \times W \times H$: 510 \times 100 \times 75 mm

Material: Silicone gel



Flat bottom chest roll medium

 Protects from positioning injuries at legs or chest

Art.-no. 61921

L x W x H: 405 x 100 x 75 mm

Material: Silicone gel



Flat bottom chest roll small

 Protects from positioning injuries at legs or chest

Art.-no. 61920

L x W x H: 305 x 100 x 75 mm

Material: Silicone gel





Bottom roll children

- For positioning of patient in supine position
- Protects from positioning injuries at legs or feet

Art.-no. 61926

L x W x H: 500 x 150 x 75 mm

Material: Artificial leather covered foam

Electrically conductive



Bottom roll adults

- For positioning of patient in supine position
- Protects from positioning injuries at legs or feet

Art.-no. 61925

L x W x H: 500 x 220 x 110 mm

Material: Artificial leather covered foam

Electrically conductive



Foot positioning pad

Protects from positioning injuries at legs

Art.-no. 61930


L x W x H: 600 x 220 x 180 mm

Material: Artificial leather covered

foam

Electrically conductive

Positioning aids



Leg positioning tunnel

- For positioning patient in lateral position
- Protects from positioning injuries at leg
- Radiolucent

Art.-no. 61935

L x W x H: 600 x 390 x 220 mm

Material: PU coated foam

Electrically conductive



Heel pads

- For positioning patient in supine position
- Protects from positioning injuries at heels
- Pair

Art.-no. 61940

L x W x H: 190 x 100 x 70 mm


Material: Silicone gel



Our full range of products:































Made in Germany